


UltraSEEK™ Oncogene Panel



Screen for Low Abundance Somatic Mutations in Heterogeneous Tumor and Blood Samples at $\leq 1\%$ Frequency

Agena Bioscience's UltraSEEK™ Oncogene Panel is a set of pre-validated assays comprising 26 driver mutations in 12 oncogenes commonly observed in neoplastic cells from various tumor types, for use on the MassARRAY® System. The panel offers a novel approach for ultrasensitive multiplexed mutation detection to aid in translational and clinical research studies.

- ✓ Compare mutation profiles for matched tissue, circulating cell-free DNA, and tumor cell DNA samples.
- ✓ Identify mutations that may not have been detected with lower sensitivity methods.
- ✓ Confirm mutation status from other molecular methods, including Agena Bioscience's OncoCarta™ Panel or next generation sequencing panels.

The UltraSEEK Oncogene Panel requires only 40 ng input DNA to assay. Only non-synonymous coding mutations previously reported to occur in human tumor samples were selected, including those in key genes implicated in resistance to tyrosine kinase inhibitors.

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@AgenaBio.com
Web: agenabioscience.com

US: +1.877.4.GENOME
EU: +49.40.899676.0
AP: +61.7.3088.1600
JP: +81.3.6231.0727
CN: +86.21.6427.0566

51-20044R4.0 0316

Agena
BIOSCIENCE






THE MASSARRAY WORKFLOW

The UltraSEEK Oncogene Panel consists of multiplex PCR followed by a mutation-specific single base extension reaction. The UltraSEEK Oncogene Panel requires 4 wells per sample tested. In addition to the mutation-specific markers, each well contains six controls to gauge the performance of the process, including bead capture, clean-up, and elution steps. The extension reaction utilizes a single mutation-specific chain terminator labeled with a moiety for solid phase capture. After capture, clean-up, and elution, products are spotted onto a SpectroCHIP® Array for analysis on the MassARRAY System.

THROUGHPUT

The UltraSEEK Oncogene Panel contains multiplexed assays in 4 wells. The panel can be run in 24-well format (6 samples per plate), 96-well format (24 samples per plate), or 384-well format (96 samples per plate). Six to several hundreds of samples can be processed per day, providing flexibility in sample throughput and batching requirements.

PANEL COMPONENTS†

AMPLIFY		PCR Enzyme PCR Accessory Set UltraSEEK™ Oncogene PCR Primers
EXTEND		UltraSEEK Accessory Set UltraSEEK Termination Mixes UltraSEEK Oncogene Extend Primers
CAPTURE		UltraSEEK Capture Reagent Set Streptavidin Beads and Buffers
DETECT		SpectroCHIP® Array and Clean Resin (1 x 384, 4 x 96, 2 x 96)
ANALYZE		MassARRAY® Analysis Software

†A plate magnetic stand is required but not included. Please contact Agena Bioscience for recommendations.

ULTRASEEK ONCOGENE PANEL GENES AND MUTATIONS

GENE	AA
<i>ABL1</i>	T3151
<i>AKT1</i>	E17K
<i>ALK</i>	F1174L
<i>BRAF</i>	V600E
	V600K
	V600R
<i>EGFR</i>	G719S
	T790M
	L861Q
	E746-A750del
<i>FLT</i>	1836del
<i>IDH1</i>	R132H
<i>IDH2</i>	R140Q
	R172K
<i>JAK2</i>	V617F
<i>KRAS</i>	G12A
	G12D
	G12V
	G12C
	G12S
<i>NRAS</i>	G13D
	Q61K
<i>PIK3CA</i>	Q61R
	E545K
	E542K
	H1047R

ORDERING INFORMATION

CAT NO	FORMAT	SAMPLES/KITS
27600	1 x 384	96
27601	4 x 96	96
27602	2 x 96	48
27603	8 x 24	48

For Research Use Only. Not for use in diagnostic procedures.

MassARRAY and SpectroCHIP are registered trademarks of Agena Bioscience, Inc. Agena Bioscience, UltraSEEK, and OncoCarta are trademarks of Agena Bioscience. ©2016 Agena Bioscience, Inc. All rights reserved.